EN Mediblink Portable mesh nebulizer Alien M425

INSTRUCTIONS FOR USE

To ensure proper operation of this product, please read this manual carefully before use.

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TABLE OF CONTENTS

1. Product Description	7
2. Safety Precautions	8
3. Structure and Composition of the Product	15
4. Names and Functions of Components	16
5. Power Supply	17
6. Assembly and Medication Filling	20
7. Nebulization Procedure	22
8. Cleaning and Disinfection	24
9. Troubleshooting	27
10. Specifications	29
11. EMC Technical Information	31
12. Product Warranty	36

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- Thank you for purchasing this product.
- To ensure proper operation of this product, please read this manual carefully before use.
- It is recommended to keep this manual in a place where it can be easily accessed at any time.
- The illustrations contained in this manual are for reference only.
- Please follow the guidelines provided in this manual for installation and use of the product.
- This is a medical device; when using this equipment, please strictly follow your doctor's recommendations.

1. Product Description

Information to
Know Before Use

Intended Use

Mediblink Portable mesh nebulizer Alien is designed for the delivery of aerosolized medications for therapeutic purposes.

Intended Operators

- Qualified medical professionals such as doctors, nurses, and therapists.
- Caregivers and patients following the guidance of medical professionals.
- Individuals capable of understanding the basic operation of the device and its user manual.

Contraindications

- Not to be used by individuals with acute cardiopulmonary failure or airway obstruction.
- Not suitable for individuals allergic to aerosolized medications.

Operating Environment

This product is suitable for use in medical facilities such as hospitals, clinics, and doctor's offices, as well as standard home environments.

Precautions for Use

- This product is a medical device and must be used under the guidance of a doctor. Ensure that the type of medication, dosage, and method of use strictly follow the doctor's instructions.
- Note that this product is not designed for general nasal or pharyngeal mucosa spray humidification.

- The nebulization performance may vary depending on the properties of the medication. Specifically, liquids with higher surfactant activity or viscosity, such as solvents and expectorants, may result in slower nebulization rates. Additionally, lower medication temperatures may also decrease the nebulization rate.
- This product is suitable for home use and can be used by both adults and children.
- Do not use purified water as the nebulizing solution.

2. Safety Precautions

Information to
Know Before Use

Before using this product, please carefully read this user manual.

- **Warnings and Symbols:**

The warnings and symbols included in this manual are intended to ensure safe and proper operation of the product, minimizing the risk of harm to you or others.

- **Explanation of Warning Symbols and Icons:**

Below are the warning symbols, icons, and their respective explanations:

About Warnings and Cautions	
 Warning	Improper use of the product may result in death or serious injury.
 Caution	Improper use may cause personal injury or damage to property.

Symbol Descriptions

	BF-Type Applied Part		Date of Manufacture
	Before use, read the instructions. Electronic instructions for use: http://www.mediblink.com/f/m425.pdf		Non-Ionizing Radiation
	Serial Number		Medical Device
	Manufacturer		MR Unsafe
	Authorized Representative in the European community		Unique Device Identifier
	Warning!		Batch number*
	Reusable for Single Patient		Product reference number
	Do Not Dispose Randomly		Model number
	CE Marking of Conformity, Notified body ID number		Ingress protection degree

*The production date can be read from LOT number [YYYYMM]; first four digits represent the Year and last two digits Month of production. Example: LOT 202503 = March 2025

**WARNING:*****a) While in Use***

This product must be purchased and used under the guidance of a physician. Follow your doctor's recommendations regarding the type of medication, dosage, and usage method.

If you experience allergic reactions or any other issues while using this product, stop using it immediately and seek medical assistance.

Do not use this product if you have acute cardiopulmonary failure or airway obstruction. Consult a physician for appropriate advice.

Before the first use, after long periods of non-use, or if components become contaminated, it is essential to clean and disinfect all parts thoroughly.

Please avoid applying any oily skin creams before using this product and ensure all oral secretions and food residues are removed prior to use. After use, promptly clean your face.

This product should only be used by children or individuals with special needs under proper guidance and supervision.

Discard any used medication solution after each use, always use fresh medication for each inhalation session.

- Reuse of medication may result in bacterial infections, potentially worsening symptoms.

This product is designed for use by a single individual. Sharing the device with others is prohibited.

Do not use water for nebulization, as it may worsen symptoms.

This product is designed exclusively for therapeutic purposes. Do not use it for any other purpose. Misuse carries risks, and the manufacturer is not responsible for any consequences arising from improper use.

Avoid contact with the nebulizer mesh. As this may damage the nebulizer mesh and potential do harm to the patient.

This product should not be used by unconscious patients or those unable to breathe independently.

Do not use or store this product in environments with flammable gases or high humidity, such as bathrooms. Doing so may result in fire or electric shock.

Avoid placing the product in environments with frequent temperature and humidity fluctuations, such as inside a car during summer or in direct sunlight.

Do not use the product in anesthesia circuits, ventilatory breathing systems, or oxygen supply pipelines. As this may lead to severe accidents.

Do not modify this device.

Replacement of components should only be performed by qualified technicians approved by the manufacturer.

Ensure the product is not placed in areas where it is difficult to access or operate the power disconnect mechanism.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

b) Power

Do not use damaged or cracked power cords or plugs.

- Using such equipment may result in electric shock, short circuits, or fire hazards.

Ensure that the power adapter used is compatible with AC 100V-240V. Do not use multiple power cords simultaneously, and each socket should only connect to one device.

- This precaution helps prevent fire, electric leakage, electric shock, or device malfunction.

Do not plug in or unplug the power cord with wet hands.

- Doing so may lead to electric shock or injury.

If the device is not used for an extended period (more than 3 months), remove the batteries. Replace depleted batteries with two new ones of the same type immediately.

- This prevents potential damage to the device caused by battery leakage.

Do not mix different types of batteries or use a combination of new and old batteries.

- Do not use lithium batteries.

Install the batteries according to the polarity markings on the device, ensuring the positive and negative ends are correctly aligned.

- Incorrect installation may cause batteries to overheat, leak, or rupture, leading to damage to the device.
- * The power adapter must comply with the requirements of the IEC 60601-1 standard.

c) Maintenance and Storage

Ensure all components are stored in a clean environment to prevent infection. Cleaning and disinfecting the device are necessary before first use, after prolonged disuse, or when components are contaminated.

- Neglecting these steps may result in bacterial infection or worsening of symptoms.

Always clean the medication cup, mouthpiece, and mask after each use. Once cleaned, disinfect and dry them immediately.

- Failing to do so may promote bacterial growth.
- Failing to do so may clog the nebulizer plate, preventing proper mist generation.
- Failing to do so may cause new medication to mix with residual medication, potentially exacerbating symptoms during inhalation.

Ensure the power adapter is unplugged before starting the cleaning process.

- This helps prevent electric shock or injury.

Do not use cotton swabs, brushes, or hard objects to clean the nebulizer plate.

Avoid spilling water or medication on the main unit or power adapter. Do not immerse them in water or other liquids. If liquid is accidentally spilled, unplug the power adapter immediately and clean it up.

- This precaution helps prevent damage, electric shock, or device malfunction.

Store the device in a location out of the reach of infants and children.

- Small components may pose a choking hazard. In such cases, seek medical attention immediately.

WARNING:

USAGE PRECAUTIONS

a) While in Use

Ensure that no mobile phones or other electromagnetic radiation sources are within at least 30 cm of this device during use.

- Such devices may interfere with performance and reduce efficiency.
- Avoid prolonged close contact between the mask and your face. This could lead to skin irritation.

This product is intended for human use only and must not be used on other animals.

Avoid allowing nebulized medication to come into contact with your eyes.

- Excessive or improper exposure to medication may cause eye irritation.
- Ensure that the medication cup contains liquid before using this product to avoid damaging the nebulizer plate.

Do not exceed the maximum capacity of the medication cup.

- Overfilling may cause leakage, increasing the risk of damage, electric shock, or device malfunction.

Clean any stains on the electrodes of the medication cup or main unit promptly.

- Dirty electrodes may reduce the product's efficiency.

Take care to prevent the main unit or its components from falling or experiencing strong impacts.

- Such incidents may result in damage, electric shock hazards, or malfunction.

During use, do not cover the device with bedclothes, blankets, or towels.

- Doing so may lead to overheating, electric shock, or device damage.

Do not wrap the power cord around your neck, as this could result in accidents or injuries.

b) Maintenance and Storage

Do not heat this product or its components in a microwave, place them in a dishwasher, use a hair dryer to dry them, or boil them.

- High temperatures may cause deformation, affecting the normal nebulization function.
- This may also pose a fire hazard.

Ensure no detergent residues remain on the components. After disinfecting with a disinfectant, rinse the components thoroughly with clean water to remove all disinfectant residues.

- Residual disinfectant may worsen symptoms when using the device.

If the device will not be used for an extended period, remove the batteries from the battery compartment.

- Used batteries must be disposed of according to local regulations.
- Do not carry or store the nebulizer with residual medication still inside.

* This product may not meet performance specifications if used or stored in temperature and humidity conditions not recommended by the manufacturer, for details on the temperature and humidity ranges for operation and storage, please refer to the specifications section.

3. Structure and Components of the Product

Information to
Know Before Use

The product consists of the following components:

Main unit, Medication cup, Adult mask, Child mask, Mouthpiece, USB-C cable.

The packaging box includes the components listed below.
If any items are missing, please contact us immediately.

Configuration list


	Main unit		Medication Cup (Composed of a quick-release nebulization nozzle and a liquid storage component.)
	1 adult mask		1 child mask
	1 mouthpiece		1 USB-C cable
	1 user manual		

4. Names and Functions of Components

Information to
Know Before Use

Medication Cup

Includes a Metal Nebulizer Plate for Nebulization.


Medication Cup

Stores the Medication to be Nebulized.

Nebulizer plate

The liquid medication is atomized through the tiny pores in the center of the high-frequency vibrating nebulizer plate, forming a fine mist.

Main unit

Electrodes

Designed to connect with the medication cup, creating a closed circuit to enable the nebulization function.

Power Button

Used to turn the nebulizer on, adjust the nebulization rate, and turn it off.

Indicator Light


Displays the current status of the nebulizer.

Medication Cup Release Button

Press to disconnect the medication cup from the main unit.

USB Port

Connects to the USB cable, allowing the main unit to be powered through this interface.


5. Power

Information to
Know Before Use

Both the AA battery version and lithium battery version of this product can be use directly by connecting the power through USB cable.

Use the product with USB cable

1. Connect one end of the USB cable to the USB interface of the host, Plug the other end into the USB power adapter.
2. Plug the USB adapter into a power outlet.

Regular AA battery version

This product uses two 1.5V AA batteries as its power source please always check the battery condition before each use.

Install or change new batteries

Open the battery cover

Gently open the battery compartment cover at the bottom.

Insert the batteries

Insert the batteries following the polarity markings inside the battery compartment.

Close the battery cover

Close the battery compartment cover. Press firmly until you hear a “click”, indicating that the cover is securely fastened.

Battery Status Indicator

Blue Light Flashing

When the battery power is low, the indicator light will flash blue. Replace both batteries with new ones as soon as possible.

- Do not short-circuit the batteries or install them with reversed polarity during installation and use, as it may cause severe damage to the device.
- If the device will not be used for an extended period, remove the batteries from the battery compartment.
- Batteries should be treated as chemical waste and disposed of in accordance with local regulations.

Lithium battery Version

1. This product is powered by, built-in rechargeable lithium battery.
2. This lithium battery is rechargeable, can be charged by connecting to the power through the USB cable. During the charging process, the indicator light will flash in color blue. When the battery is fully charged, the indicator light will be stable in color blue.
3. Lithium battery can work for at least 240 minutes after fully charged.
4. The capacity of the Lithium battery will still be greater than 80% after re-charges over 300 times.

NOTE:

1. Before nebulization, check if the battery is sufficiently charged.
2. When the low battery indicator lights up, recharge the battery promptly.
3. To ensure the normal life span of the lithium battery, it should be charged at least once a month; if not used for an extended period, fully charge the battery before reuse.
4. Please dispose of the used battery according to local regulations.
5. To charge the lithium battery, please use the power adapter that meets the requirements of the IEC60601-1 standard.

6. Assembly of Components and Adding Medication

Information to
Know Before Use

WARNING:

Before the first use, after prolonged disuse, or if components are contaminated, ensure that all components are thoroughly cleaned and disinfected.

Before Use

1. Carefully read the instruction manual.
2. Disassemble, clean, and disinfect all components.
3. Assemble the components only after they are cleaned, disinfected, and dried.

Adding Medication

Open the medication cup lid

- Push the latch at the back of the medication cup lid outward to unlock it.
- Lift the medication cup lid upward.

Tip: When unlocking the latch, gently press down on the medication cup lid to make the operation easier.

Add medication

- If the medication requires dilution or mixing, prepare it in a clean, disinfected container beforehand.
- Pour between 0.5 mL and 8 mL of medication into the medication cup. Ensure the medication temperature does not exceed 60°C.

Close the medication cup lid

- Lower the medication cup lid back into place.
- Securely latch the back of the lid. A “click” sound indicates that it is locked.

Tip: When locking the latch, gently press down on the medication cup lid to make the operation easier.

Note: Avoid using thick, high-viscosity, or highly concentrated liquids, as these may clog or damage the nebulizer components, affecting normal functionality.

Assembly of Components

Installing the Medication Cup

- Press and slide the medication cup onto the main unit in the direction indicated in the diagram.
- When you hear a “click”, it means the medication cup is securely installed.


Installing the Mask or Mouthpiece

1. Installing the Mouthpiece

- Choose either the adult mask or child mask based on the user's needs.
- Attach the selected mask to the nebulization outlet of the medication cup.

2. Installing the mask

- Do not touch the nebulizer plate on the medication cup.
- Use the adult mask for adults and the child mask for children.
- Ensure that only masks and mouthpieces specified by the manufacturer are used.
- When disassembling, perform the steps in reverse order.

7. How to Perform Nebulization

Information to
Know Before Use


WARNING:

Follow your physician's recommendations for the type of medication, dosage, and method of use.

1. Adjusting the Main Unit Position


**Upright
Position**

(✓)


**Tilting
Downward**

(✓)

**Tilting
Upward**

(✗)

Before starting nebulization, ensure the main unit is either upright or tilted downward. This ensures continuous contact between the medication and the nebulizer plate.

2. Starting Nebulization


1. Press the power button briefly.

The indicator light will flash green, and the nebulizer enters the auto-cleaning mode.

This mode intermittently sprays mist to rinse the nebulizer plate.

2. After auto-cleaning, the indicator light will remain solid green, indicating the first-level nebulization speed.

- While inhaling, remain calm and relaxed. Take slow, deep breaths to ensure the medication reaches deep into your lungs.
- Hold your breath briefly after inhaling, and exhale slowly once the mouthpiece or mask is removed.
- Avoid breathing too quickly. If necessary, take short breaks between nebulizations.

Reminder:

- High-concentration medication may reduce the nebulization speed.
 - If excess medication accumulates on the nebulizer plate, nebulization may be interrupted. In this case:
 - Turn off the power.
 - Remove the mouthpiece or mask.
 - Use a clean cloth to absorb the excess liquid on the nebulizer plate.
3. To increase the nebulization speed, briefly press the power button again. The indicator light will turn solid blue, indicating the second-level mode.
 4. To further increase the nebulization speed, briefly press the power button again. The indicator light will turn solid cyan, indicating the third-level mode.

3. End the Nebulization

In the third-level nebulization mode, press the power button briefly again to stop nebulization.

Reminder:

- The nebulizer's single session duration is 10 minutes. After 10 minutes, the device will automatically shut off. To continue using it, restart the device.
- When the medication in the cup is fully nebulized, the nebulizer will activate its waterless detection feature. The indicator light will flash yellow to signal this, after which the device will stop working and power off automatically.

8. Cleaning and Disinfection

Information to
Know Before Use


WARNING:

After each use, ensure the medication cup, mouthpiece, and mask are thoroughly cleaned. Once cleaned, immediately disinfect and air-dry them. Fail to clean and disinfect components properly and regularly may result in bacterial growth, increasing the risk of infection.

Note:

When cleaning and disinfecting the medication cup, do not disassemble the quick-release nebulization nozzle. Improper disinfection or cleaning may cause rust or damage to the components. The quick-release nebulization nozzle is designed for replacement only.

Removing the Medication Cup

1. Ensure the main unit is turned off, the USB cable is disconnected, and the USB silicone cover is in place.
2. Press the medication cup release button on the back of the main unit while pulling the medication cup upward to remove it.

Disposing of Residual Medication

Empty any remaining medication from the medication cup.


Nebulization Cleaning with Water

- Refill the medication cup with water (do not use purified water) and close the lid.
- Reattach the medication cup to the main unit and nebulize for 1-2 minutes.
- After nebulization, remove the medication cup and pour out the remaining water.

Reminder:

- This step effectively cleans the nebulizer plate and extends its lifespan.
- Do not inhale the water during nebulization.

Cleaning and Disinfecting the Medication Cup, Mouthpiece, and Mask

- Fill a clean container with purified water at a temperature not exceeding 40°C.
- Open the medication cup lid.
- Immerse the medication cup, mouthpiece, and mask in the water.
- After 20 minutes of soaking, remove the components and wipe them clean with sterile medical gauze.
- Replace the purified water with 75% medical-grade alcohol and repeat the above steps for disinfection.

Cleaning and Disinfecting the Medication Cup, Mouthpiece, and Mask

- Wipe the surface of the main unit with a damp piece of medical-grade gauze.
- Wet another piece of gauze with 75% medical-grade alcohol and wipe the main unit surface.
- Use dry medical-grade gauze to wipe the unit until no moisture remains on the surface.

Disassembly and Assembly

Replacing the Quick-Release Nebulization Nozzle.

Disassembly

- Rotate the nozzle counterclockwise by 25 degrees until the alignment marks on the nozzle and medication cup are aligned.
- Pull the nozzle out in the direction indicated by the arrow.

Assembly

- Align the marks on the nozzle and medication cup, and insert the nozzle into the medication cup following the direction of the arrow.
- Rotate the nozzle clockwise by 25 degrees until you hear a “click”, indicating it is securely in place.

9. Troubleshooting

If any
Questions

When encountering an issue, first ensure there are no electronic devices within 30 cm that could cause interference.

If the problem persists, refer to the table below.

Issue	Possible Cause	Solution
Nebulization rate is very low.	Low battery (indicator light flashing blue).	Power supply the device via the USB port or replace the batteries.
	Nebulizer plate is dirty or clogged.	Clean and disinfect the nebulizer plate. Replace the medication cup.
	Insufficient medication volume.	Ensure enough medication is added, but do not exceed 8 mL.
Nebulization rate is very low.	Medication does not meet nebulization requirements.	Use appropriate medication that is water-soluble and non-corrosive. If dilution is needed, follow your physician's advice before adding it to the medication cup.
	Water detection probe in the medication cup is dirty.	Thoroughly clean and disinfect the medication cup to restore the probe to its original clean state.
	Incorrect handling angle of the device.	Ensure the device remains upright or slightly tilted downward so the medication is in full contact with the nebulizer plate.

The nebulizer does not operate when turned on.	Batteries are inserted incorrectly, or battery power is low.	Re-insert the batteries correctly or replace them.
The nebulizer does not operate when turned on, and the indicator light flashes blue/yellow.	Low battery.	Power supply the device via the USB port or replace the batteries.
	Insufficient medication volume, triggering the automatic waterless detection shutdown.	Add enough medication, but do not exceed 8 ml.
	Medication cup is not securely attached to the main unit.	Reinstall the medication cup onto the main unit until you hear a "click."
The nebulizer operates but shuts down automatically after a while.	Low battery (indicator light flashing blue).	Power supply the device via the USB port or replace the batteries.
	Medication cup is not securely attached to the main unit.	Reinstall the medication cup onto the main unit until you hear a "click."
	Insufficient medication volume, triggering the automatic waterless detection shutdown.	Add enough medication, but do not exceed 8 ml.
	The nebulization session reaches 10 minutes, triggering the automatic shutdown.	This is normal. If continued nebulization is required, press the power button again.
	Medication is not in full contact with the nebulizer plate.	Gently shake the nebulizer to ensure proper contact between the medication and the nebulizer plate.

Note: If the above causes are checked and your device still does not function properly, please contact distributor for assistance.

10. Specifications

Product Information

Product Name	Nebulizer (Ultrasonic mesh type)
Model/Reference number	LT-N100/M425
Power Supply	2 X 1.5 V AA batteries (AA battery version), built-in 3.7 V batteries (Lithium battery Version) or DC 5 V, 1 A (Class II adaptor with IEC 60601-1 approval)
Power Consumption	< 3.0 W
Ultrasonic Vibration Frequency	110 KHz $\pm$ 10 %
Maximum Nebulization Rate	$\geq$ 0.2 ml/min
Medication Cup Temperature	$\leq$ 60°C
Noise Level	$\leq$ 50 dB (A-weighted)
Median Particle Diameter (MMAD)	2.2 μ m $\pm$ 25%
Volume Distribution of Equivalent Particle Diameter	Proportion of particles with diameters of 1 μ m to 5 μ m $\geq$ 78%
Operating Mode	Continuous operation
Continuous Working Time	10 minutes for each use, 3 times consecutive nebulization
Maximum Medication Cup Capacity	8 ml
Dimensions	Approx. Length 74 mm $\times$ Width 50 mm $\times$ Height 127 mm (without mask or mouthpiece)

Weight	Approx. 100 g (AA battery version) Approx. 140 g (Lithium battery Version)
Ingress Protection Rating	IP22
Operating Environment	Temperature: 5°C to 35°C Relative Humidity: ≤ 80%, no condensation Atmospheric Pressure: 86.0 kPa to 106.0 kPa
Transport/Storage Environment	Temperature: -20°C to 55°C Relative Humidity: ≤ 93%, no condensation Atmospheric Pressure: 70.0 kPa to 106.0 kPa
Electromagnetic Compatibility Classification	Group I, Class B
Electric Shock Degree Classification	BF-type applied part
Electric Shock Type Classification	Class II device, internally powered device
Safety Classification	Equipment unsuitable for use in the presence of flammable anesthetic gas mixed with air, oxygen, or nitrous oxide
Disinfection/Sterilization Methods	Follow the manufacturer's recommended procedures

Note:

- The product may not operate properly if temperature or voltage conditions differ from those stated in the manual.
- Median particle diameter and equivalent particle diameter distribution were measured under the following conditions: 0.9% sodium chloride solution, $23 \pm 2^\circ\text{C}$, and relative humidity of 56%-59%. Results may vary based on medication and environmental conditions.

11. EMC Technical Information

Product
Information

Note:

1. This instrument meets the requirements of IEC60601-1-2 standard for electromagnetic compatibility.
2. The user shall install and use the device according to the EMC information provided in this file.
3. Portable and mobile RF communication equipment may affect the performance of the instrument, avoid strong electromagnetic interference when using, such as close to the mobile phone, microwave oven, etc.
4. The guidance and manufacturer's declaration are detailed in the table below.

WARNING:

1. The M425/LT-N100 medical ultrasonic nebulizer should not be used close to or stacked with other equipment. If it must be used near or stacked, it's important to observe and verify that it operates normally under the specific configuration being used.
2. Apart from the cables provided by the instrument manufacturers as replacement parts for internal components, using alternative accessories and cables might lead to higher emissions or decreased immunity.

Cable Information	Cable Length (m)	Shielded
USB Cable	0.7 m	NO

Guidance and Manufacturer Declaration – Electromagnetic Emissions

The instrument is intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in such electromagnetic.

Emission Test	Conformity	Electromagnetic environment – guidance
Radio frequency emission EN 55011	Group 1	The instrument uses RF energy only for its internal functions. Therefore, its RF emissions are low and there is little possibility of interference with nearby electronic equipment.
Radio frequency emission EN 55011	Class B	The instrument is suitable for use in all facilities, including domestic facilities and direct connection to residential low-voltage power grids.
Harmonic emission IEC 61000-3-2	N/A	
Voltage fluctuation / flicker emission IEC 61000-3-3	N/A	

Guidance and Manufacturer Declaration – Electromagnetic Immunity

The instrument is intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in such electromagnetic environment:

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge IEC 61000-4-2	±8 kV contact discharge ±2 kV, ±4 kV, ±8 kV, ±15 kV air discharge	±8 kV contact discharge ±2 kV, ±4 kV, ±8 kV, ±15 kV air discharge	The ground should be wood, concrete or ceramic. If the floor is covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient burst IEC 61000-4-4	±2 kV to power line, 100 KHz frequency	N/A	/
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to ground	N/A	/
Voltage dip IEC 61000-4-11	0 % U_T ; lasts 0.5 cycle 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U_T ; 1 cycle 70 % U_T ; 25/30 cycles	N/A	/
Voltage interruption IEC 61000-4-11	0 % U_T ; 250/300 cycles	N/A	/
Power frequency magnetic field IEC 61000-4-8	30 A/m	30 A/m, 50/60 Hz	/

Note: U_T refers to the AC grid voltage before the test voltage is applied.

Guidance and Manufacturer Declaration – Electromagnetic Immunity

The instrument is intended to be used in the electromagnetic environment specified below, and the purchaser or user should ensure that it is used in such electromagnetic.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Radio frequency conduction IEC 61000-4-6	3 V 150 kHz ~ 80 MHz 6 V 150 kHz ~ 80 MHz 80% AM, 1 kHz	3 V 6 V	Portable and mobile RF communications equipment should not be used closer to any part of the instrument than the recommended isolation distance, including cables. This distance is calculated by a formula corresponding to the transmitter frequency. Recommended isolation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz ~ 800 MHz $d = 1.2 \sqrt{P}$ 800 MHz ~ 2.7 GHz
Radio frequency radiation IEC 61000-4-3	10 V/m 80 MHz ~ 2.7 GHz 80% AM 1 kHz	10 V/m	In the formula: P – Maximum rated output power according to the transmitter manufacturer, in watts (W); d – The recommended isolation distance, in meters (m). The field strength of a fixed RF transmitter is determined by a survey c of the electromagnetic field, and each frequency range d should be lower than the compliance level. Interference may occur near devices marked with the following symbol. 

Note 1: At 80 MHz and 800 MHz frequencies, the formula for the higher frequency band is used.

Note 2: These guidelines may not be suitable for all situations. Electromagnetic propagation is affected by the absorption and reflection of buildings, objects, and human bodies.

- A For fixed transmitters, such as base stations for wireless (cellular/cordless) telephones and terrestrial mobile radios, amateur radios, AM and FM radio broadcasts, and television broadcasts, the field strengths are not predictable theoretically. To assess the electromagnetic environment of a fixed RF transmitter, the survey of the electromagnetic field should be considered. If the field strength of the location where the instrument is located is higher than the above-mentioned applicable RF compliance level, the instrument should be observed to verify that it is functioning properly. Additional measures may be necessary if abnormal performance is observed, such as reorienting the instrument.
- B The field strength should be less than 3V/m over the entire frequency range from 150 kHz to 80 MHz.

Recommended Isolation Distance Between Portable and Mobile RF Communications Equipment and Instruments

The instrument is expected to be used in an electromagnetic environment where radio frequency disturbances are controlled. Depending on the maximum rated output power of the communication device, the purchaser or user can prevent electromagnetic interference by maintaining the minimum distance between the portable and mobile RF communication device (transmitter) and the instrument as recommended below.

Rated maximum output power of the transmitter/ W	Isolation distance for different frequencies of the transmitter/m		
	150 KHz ~ 80 MHz $d = 1.2 \sqrt{P}$	80 MHz ~ 800 MHz $d = 1.2 \sqrt{P}$	800 MHz ~ 2.7 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For the maximum rated output power of the transmitter not listed in the above table, the recommended isolation distance d (in meters (m)) can be determined by the formula in the corresponding transmitter frequency column, where P is the maximum rated output power provided by the transmitter manufacturer, in watts (W).

Note 1: At the 80 MHz and 800 MHz frequency points, the formula for the higher frequency band is used.

Note 2: These guidelines may not be suitable for all situations. Electromagnetic propagation is affected by the absorption and reflection of buildings, objects and the human body.

12. Product Warranty

Product
Information

Product: Mediblink Portable mesh nebulizer Alien M425

Manufactured for (Importer for EU and Distributor): Mediblink d.o.o.,
Gubčeva cesta 19, 8210 Trebnje, Slovenia; info@mediblink.com;
www.mediblink.com

Seller's name, address, signature and stamp*:

Date of extradition/sales*:

*If the invoice is accompanied by this warranty, and if all above information can be seen from the invoice, it is not necessary to fill in this field.

Warranty terms

Dear customers!

The warranty period is **3 years** and starts from the day of product purchase. The warranty period on spare parts: Adult mask, Child mask, Medication cup, Mouth-piece, Mesh membrane is **1 year**. In case of product claim, you have to show the invoice. We kindly ask you to save the invoice!

Unfortunately, wrong handling with the device is a reason for 95% of customer complains. You can easily avoid any problem, by getting useful information provided by our special service department. To reach our service department, you can call or send e-mail to Mediblink local Distributor.

Before sending the product back to retailer, we kindly ask you to call our service department, to get help about how to use the device to save you with unneeded trips.

The Manufacturer guarantees free elimination of all imperfections due to defects in material or manufacturing procedure by repairing or replacing the product. In case that the product cannot be repaired

or replaced, the customer will get the money refund. The guarantee is not valid in case of the force majeure, accidents or unexpected events (such as lightning, water, fire etc.), incorrect use or incorrect transport, non-compliance with safety and maintaining regulations or in case of unprofessional product intervention.

Traces of every day product usage (scratches, abrasions) and not subject to claim. The warranty does not eliminate the customer rights, which originate from seller responsibility for product flaws. By accepting the claimed product by the service department, the service department does not take responsibility for loss of saved data or settings on the product. All product repairs, which are performed out of product warranty period, have to be paid by customer by prior notice.

The Manufacturer guarantees the product quality and flawless product operation in the warranty period, which starts with the day of product purchase. If the product cannot be repaired in 45 days, the product will be replaced with a new one. In case that the product cannot be replaced, the money will be refunded to the customer.

In case of product claim, call or send e-mail to Mediblink local distributor.

Replacement parts can be ordered at www.mediblink.com.

Manufacturer:
Zhuhai Linte Medical Instrument Co., Ltd.
501, Building #3, 2 Jinliang Road,
Hongqi Town, Jinwan District, 519090
Zhuhai, Guangdong,
P. R. China
E-mail: lintemed@lintemed.com

Manufactured for
(importer for EU &
Distributor):
Mediblink d.o.o.,
Gubčeva cesta 19,
8210 Trebnje, Slovenia
info@mediblink.com;
www.mediblink.com

Shanghai International Holding Corp.
GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
E-mail: shholding@hotmail.com

Instructions for
use version: A/4
Issue date: 1.10.2018
Date of last change:
15.7.2025

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